


# REUSABLE LARYNGOSCOPE SYSTEMS



Timesco: The U.K.'s No.1  
Laryngoscope Brand

**TIMESCO**  
INSTRUMENTS FOR LIFE

Sirius-XL Eclipse blades feature a flexible tip designed for difficult intubations and are available in Macintosh 4, 3 and 2. Supplied as a blade only in a case or as a blade and handle set in a case.



The non-replaceable glass fibre bundle is fully enclosed in the blade which allows easy cleaning.

Rigid metal construction from fully autoclaveable, antimagnetic stainless steel.

Lever positioned at the rear of the blade to control the tilting tip.

Fully autoclaveable metal blade, covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.

Clear product information etched on the blade.

Atraumatic distal end.

Tilting tip for elevation of the epiglottis, providing an improved view of the larynx, without excessive force.

Green ISO 7376 fitting metal block compatible with reusable handles (pages 6-9).

3 metal ball bearings for secure fitting.



CE marked



Fully autoclaveable



Latex Free



Extensive guarantees



ISO 7376 compliant



Supplied with multilingual instructions



Compatible fibre optic laryngoscope handles are colour coded

**A** | 2955.175.25 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.4

- Large Adult (No.4)
- 160 (L) x 24mm (H)

**B** | 2955.175.20 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.3

- Adult (No.3)
- 135 (L) x 24mm (H)

**C** | 2955.175.15 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.2

- Child (No.2)
- 115 (L) x 24mm (H)



## Sets in hard cases

Also available as sets in hard cases:

- 2955.176.20 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.3 WITH XENON HANDLE
- 2955.176.25 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.4 WITH XENON HANDLE
- 2955.177.25 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.3 AND NO.4 WITH XENON HANDLE

- 2955.178.20 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.3 WITH XLED HANDLE
- 2955.178.25 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.4 WITH XLED HANDLE
- 2955.179.25 | SIRIUS-XL ECLIPSE MACINTOSH BLADE NO.3 AND NO.4 WITH XLED HANDLE